

Compound TN12639

Chemical Properties


CAS No. :

Formula: C39H62O13

Molecular Weight: 738.91

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(1R,5aS,7S,9aR,11aR)-7-[[[(2R,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-[[[(2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy]oxan-2-yl]oxy]-1-[(1S)-1-[(2R,5S,6S)-6-hydroxy-5-methyloxan-2-yl]ethyl]-9a,11a-dimethyl-1H,2H,3H,3aH,5H,5aH,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-5-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3533 mL	6.7667 mL	13.5334 mL
5 mM	0.2707 mL	1.3533 mL	2.7067 mL
10 mM	0.1353 mL	0.6767 mL	1.3533 mL
50 mM	0.0271 mL	0.1353 mL	0.2707 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481