

Compound TN12644

## Chemical Properties


CAS No. :

Formula: C39H62O14

Molecular Weight: 754.91

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(1R,7S,9aR,11aR)-7-[[[(2R,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-[[[(2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy]oxan-2-yl]oxy]-3a-hydroxy-1-[(1S)-1-[(2R,5S)-6-hydroxy-5-methyloxan-2-yl]ethyl]-9a,11a-dimethyl-1H,2H,3H,3aH,5H,5aH,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-5-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3247 mL	6.6233 mL	13.2466 mL
5 mM	0.2649 mL	1.3247 mL	2.6493 mL
10 mM	0.1325 mL	0.6623 mL	1.3247 mL
50 mM	0.0265 mL	0.1325 mL	0.2649 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481