

## 4-[3-ethenyl-2-hydroxy-3-methyl-6-(prop-1-en-2-yl)cyclohexyl]phenol

## Chemical Properties


CAS No. :

Formula: C<sub>18</sub>H<sub>24</sub>O<sub>2</sub>

Molecular Weight: 272.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-[3-ethenyl-2-hydroxy-3-methyl-6-(prop-1-en-2-yl)cyclohexyl]phenol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6712 mL	18.356 mL	36.7121 mL
5 mM	0.7342 mL	3.6712 mL	7.3424 mL
10 mM	0.3671 mL	1.8356 mL	3.6712 mL
50 mM	0.0734 mL	0.3671 mL	0.7342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481