

## Compound TN12660

### Chemical Properties


CAS No. :

Formula: C<sub>30</sub>H<sub>46</sub>O<sub>4</sub>

Molecular Weight: 470.69

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	(2Z)-6-[(3aS,9aS,11aS)-1-hydroxy-3a,6,6,9a,11a-pentamethyl-7-oxo-1H,2H,3H,3aH,4H,5H,5aH,6H,7H,8H,9H,9aH,10H,11H,11aH-cyclopenta[a]phenanthren-1-yl]-2-methylhept-2-enoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1245 mL	10.6227 mL	21.2454 mL
5 mM	0.4249 mL	2.1245 mL	4.2491 mL
10 mM	0.2125 mL	1.0623 mL	2.1245 mL
50 mM	0.0425 mL	0.2125 mL	0.4249 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481