

Compound TN12688

Chemical Properties


CAS No. :

Formula: C₂₉H₄₀O₇

Molecular Weight: 500.63

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	6-[(1R)-1-[(5R,9aR,11S,11aS)-5,11-dihydroxy-9a,11a-dimethyl-9-oxo-1H,2H,3H,3aH,3bH,4H,5H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-1-yl]-2-hydroxyethyl]-3-(methoxymethyl)-4-methyl-5,6-dihydro-2H-pyran-2-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9975 mL	9.9874 mL	19.9748 mL
5 mM	0.3995 mL	1.9975 mL	3.995 mL
10 mM	0.1997 mL	0.9987 mL	1.9975 mL
50 mM	0.0399 mL	0.1997 mL	0.3995 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481