

Compound TN12689

Chemical Properties


CAS No. :

Formula: C29H40O6

Molecular Weight: 484.63

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

6-[[[(1R)-1-[(5R,9aR,11aS)-5-hydroxy-9a,11a-dimethyl-9-oxo-1H,2H,3H,3aH,3bH,4H,5H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-1-yl]-2-hydroxyethyl]-3-(methoxymethyl)-4-methyl-5,6-dihydro-2H-pyran-2-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0634 mL	10.3171 mL	20.6343 mL
5 mM	0.4127 mL	2.0634 mL	4.1269 mL
10 mM	0.2063 mL	1.0317 mL	2.0634 mL
50 mM	0.0413 mL	0.2063 mL	0.4127 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481