

## Compound TN12709

## Chemical Properties


CAS No. :

Formula: C<sub>19</sub>H<sub>24</sub>O<sub>9</sub>

Molecular Weight: 396.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,4S,5S,6R,8S,10R)-2-(hydroxymethyl)-8,10-dimethoxy-3,9-dioxatricyclo[4.4.0.0 <sup>2,4</sup> ]decan-5-yl 4-hydroxy-3-methoxybenzoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5228 mL	12.6138 mL	25.2277 mL
5 mM	0.5046 mL	2.5228 mL	5.0455 mL
10 mM	0.2523 mL	1.2614 mL	2.5228 mL
50 mM	0.0505 mL	0.2523 mL	0.5046 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481