

## Compound TN12753

## Chemical Properties


CAS No. :

Formula: C<sub>35</sub>H<sub>48</sub>O<sub>11</sub>

Molecular Weight: 644.76

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(1S,3aR,3bS,7S,9aR,9bR,11R,11a5)-1-acetyl-1,3a,3b-trihydroxy-7-(((2R,4S,5R,6R)-5-hydroxy-4-methoxy-6-methyloxan-2-yl)oxy)-9a,11a-dimethyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-11-yl 4-hydroxybenzoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.551 mL	7.7548 mL	15.5096 mL
5 mM	0.3102 mL	1.551 mL	3.1019 mL
10 mM	0.1551 mL	0.7755 mL	1.551 mL
50 mM	0.031 mL	0.1551 mL	0.3102 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481