

Compound 1084954-73-3

## Chemical Properties


CAS No. : 1084954-73-3

Formula: C42H66O13

Molecular Weight: 778.98

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(1S,3aR,3bS,7S,9aR,9bR,11R,11a5)-1-acetyl-1,3a,3b-trihydroxy-7-(((2R,4S,5R,6R)-5-(((2S,4R,5R,6R)-5-hydroxy-4-methoxy-6-methyloxan-2-yl]oxy)-4-methoxy-6-methyloxan-2-yl]oxy)-9a,11a-dimethyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-11-yl (2E)-3,4-dimethylpent-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2837 mL	6.4187 mL	12.8373 mL
5 mM	0.2567 mL	1.2837 mL	2.5675 mL
10 mM	0.1284 mL	0.6419 mL	1.2837 mL
50 mM	0.0257 mL	0.1284 mL	0.2567 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481