

## Compound TN12762

## Chemical Properties


CAS No. :

Formula: C48H76O16

Molecular Weight: 909.12

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(1S,3aR,3bS,7S,9aR,9bR,11R,11a5)-1-acetyl-1,3a,3b-trihydroxy-7-(((2R,4S,5R,6R)-5-(((2S,4S,5S,6R)-4-hydroxy-5-(((2S,4R,5R,6R)-5-hydroxy-4-methoxy-6-methyloxan-2-yl)oxy)-6-methyloxan-2-yl)oxy)-4-methoxy-6-methyloxan-2-yl)oxy)-9a,11a-dimethyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-11-yl (2E)-3,4-dimethylpent-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.100 mL	5.4998 mL	10.9996 mL
5 mM	0.220 mL	1.100 mL	2.1999 mL
10 mM	0.110 mL	0.550 mL	1.100 mL
50 mM	0.022 mL	0.110 mL	0.220 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481