

4-methyl-13-methylidenetetracyclo[10.2.1.0<sup>1,9</sup>.0<sup>3,8</sup>]pentadeca-tetraene-2-carboxylic acid

## Chemical Properties


CAS No. :

Formula: C<sub>18</sub>H<sub>18</sub>O<sub>2</sub>

Molecular Weight: 266.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-methyl-13-methylidenetetracyclo[10.2.1.0 <sup>1,9</sup> .0 <sup>3,8</sup> ]pentadeca-3(8),4,6,9-tetraene-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7546 mL	18.773 mL	37.546 mL
5 mM	0.7509 mL	3.7546 mL	7.5092 mL
10 mM	0.3755 mL	1.8773 mL	3.7546 mL
50 mM	0.0751 mL	0.3755 mL	0.7509 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481