

## Compound TN12883

## Chemical Properties


CAS No. :

Formula: C34H62O13

Molecular Weight: 678.86

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-[(3-(((2R,3R,4R,5R,6S)-4,5-dihydroxy-6-methyl-3-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)oxan-2-yl]oxy)decanoyl)oxy]dodecanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4731 mL	7.3653 mL	14.7306 mL
5 mM	0.2946 mL	1.4731 mL	2.9461 mL
10 mM	0.1473 mL	0.7365 mL	1.4731 mL
50 mM	0.0295 mL	0.1473 mL	0.2946 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481