

Compound TN12896

Chemical Properties


CAS No. :

Formula: C₂₈H₅₀O₁₃

Molecular Weight: 594.7

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

3-[(3-[(2R,3R,4R,5R,6S)-4,5-dihydroxy-6-methyl-3-[(2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy]oxan-2-yl]oxy]octanoyl)oxy]octanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6815 mL	8.4076 mL	16.8152 mL
5 mM	0.3363 mL	1.6815 mL	3.363 mL
10 mM	0.1682 mL	0.8408 mL	1.6815 mL
50 mM	0.0336 mL	0.1682 mL	0.3363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481