

Compound TN12998

Chemical Properties


CAS No. :

Formula: C47H88O15

Molecular Weight: 893.21

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

1-(hexadecanoyloxy)-3-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy)propan-2-yl hexadecanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1196 mL	5.5978 mL	11.1956 mL
5 mM	0.2239 mL	1.1196 mL	2.2391 mL
10 mM	0.112 mL	0.5598 mL	1.1196 mL
50 mM	0.0224 mL	0.112 mL	0.2239 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481