

## Compound TN13017

## Chemical Properties


CAS No. :

Formula: C<sub>29</sub>H<sub>36</sub>N<sub>4</sub>O<sub>8</sub>

Molecular Weight: 568.63

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

5-carbamimidamido-2-{6-[(1Z)-2-carboxyeth-1-en-1-yl]-9a-methyl-3-octanoyl-2,9-dioxo-2H,7H,9H,9aH-furo[3,2-g]isoquinolin-7-yl}pentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7586 mL	8.7931 mL	17.5861 mL
5 mM	0.3517 mL	1.7586 mL	3.5172 mL
10 mM	0.1759 mL	0.8793 mL	1.7586 mL
50 mM	0.0352 mL	0.1759 mL	0.3517 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481