

Compound TN13027

Chemical Properties


CAS No. :

Formula: C₂₅H₃₀O₇

Molecular Weight: 442.51

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-hydroxy-2',2'',2''',3',5,5',8'a-heptamethyl-2'',4',6'',7',8',8'a-hexahydrodispiro[2,6-dioxabicyclo[2.2.1]heptane-7,1'-naphthalene-6',3''-pyran]-3,4',6''-trione is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2598 mL	11.2992 mL	22.5984 mL
5 mM	0.452 mL	2.2598 mL	4.5197 mL
10 mM	0.226 mL	1.1299 mL	2.2598 mL
50 mM	0.0452 mL	0.226 mL	0.452 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481