


(R)-TCO-OH

Chemical Properties

CAS No. :	39637-78-0
Formula:	C ₈ H ₁₄ O
Molecular Weight:	126.2
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	(R)-TCO-OH is a alkyl chain-based linker for PROTACs which joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one targeting an E3 ubiquitin ligase and the other binding the target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade these target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.9239 mL	39.6197 mL	79.2393 mL
5 mM	1.5848 mL	7.9239 mL	15.8479 mL
10 mM	0.7924 mL	3.962 mL	7.9239 mL
50 mM	0.1585 mL	0.7924 mL	1.5848 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481