

Br-PEG3-C2-Boc

Chemical Properties

CAS No. : 782475-37-0


Formula: C13H25BrO5

Molecular Weight: 341.24

Storage: Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Br-PEG3-C2-Boc is a PEG-based linker for PROTACs, facilitating the conjugation of two essential ligands crucial for PROTAC molecule formation, and enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, and the other targets the desired protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9305 mL	14.6524 mL	29.3049 mL
5 mM	0.5861 mL	2.9305 mL	5.861 mL
10 mM	0.293 mL	1.4652 mL	2.9305 mL
50 mM	0.0586 mL	0.293 mL	0.5861 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481