

Compound T128218

Chemical Properties


CAS No. :

Formula: C₆₄H₁₀₂O₃₀

Molecular Weight: 1351.49

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3S,4S,5R,6S)-6-[[[(2S,3R,4S,5S)-2-[(4aS,6aS,6bR,9R,10R,11S,12aR,14bS)-10-[[[(2R,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-[[[(2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy]oxan-2-yl]oxy]-11-hydroxy-9-(hydroxymethyl)-2,2,6a,6b,9,12a-hexamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricene-4a-carbonyloxy]-4,5-dihydroxyoxan-3-yl]oxy]-5-hydroxy-2-methyl-4-[[[(2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxy]oxan-3-yl] 1-methyl 3-hydroxy-3-methylpentanedioate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7399 mL	3.6996 mL	7.3992 mL
5 mM	0.148 mL	0.7399 mL	1.4798 mL
10 mM	0.074 mL	0.370 mL	0.7399 mL
50 mM	0.0148 mL	0.074 mL	0.148 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481