


WAY-605471

Chemical Properties

CAS No. : 378223-81-5
 Formula: C₁₄H₁₅NO₃S₂
 Molecular Weight: 309.4
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-605471 is a AHR antagonist with IC ₅₀ of 0.5 - 1 μM.
Targets(IC ₅₀)	AhR, Aryl Hydrocarbon Receptor

Solubility Information

Solubility	DMSO: 55 mg/mL (177.76 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2321 mL	16.1603 mL	32.3206 mL
5 mM	0.6464 mL	3.2321 mL	6.4641 mL
10 mM	0.3232 mL	1.616 mL	3.2321 mL
50 mM	0.0646 mL	0.3232 mL	0.6464 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Current Patent Assignee: MAGENTA THERAPEUTICS INC - WO2020/51207, 2020, A2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481