


DX2-201

Chemical Properties

CAS No. : 2749554-00-3
 Formula: C₁₈H₂₈N₂O₆S₂
 Molecular Weight: 432.55
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	DX2-201 is a potent and selective inhibitor of oxidative phosphorylation (OXPHOS) complex I, with an IC ₅₀ of 312 nM, and it exhibits anticancer activity.
Targets(IC ₅₀)	Others,OXPHOS,Mitochondrial Metabolism

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3119 mL	11.5594 mL	23.1187 mL
5 mM	0.4624 mL	2.3119 mL	4.6237 mL
10 mM	0.2312 mL	1.1559 mL	2.3119 mL
50 mM	0.0462 mL	0.2312 mL	0.4624 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481