

AKT2 Protein, Human, Recombinant (His)

General Information

Protein Construction:	Recombinant full length human AKT2 was expressed by baculovirus in Sf9 insect cells using a N-terminal His tag.
Species:	Human
Expression Host:	Baculovirus-Insect Cells
Accession:	P31751
Molecular Weight:	~60 kDa

QC Testing

Representative specific activity data for AKT2 were obtained using different assay formats:

Biological Activity:	- 44 nmol/min/mg, determined by radioactive kinase assay - 29-40 nmol/min/mg, determined by ADP-Glo kinase assay
----------------------	---

Different lots may be tested using different assay methods; therefore, values obtained using different assay formats are not directly comparable.

Purity: >70% as determined by SDS-PAGE.

Formulation: Supplied as sterile 50 mM sodium phosphate, pH 7.0, 300 mM NaCl, 150 mM imidazole, 0.1 mM PMSF, 2 mM DTT, 25% glycerol.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Enzymes are highly recommended to be shipped at frozen temperature with dry ice. Shipment made at ambient temperature may seriously affect the activity of the ordered products.

Protein Background

AKT (AK mouse plus Transforming or Thymoma) is a frequent oncogene expressed in most tissues which includes three isoforms AKT1, AKT2, and AKT3. Hyperactivation of AKT signaling is a central key in many human cancer progressions, through modulating angiogenesis, tumor growth, and cell migration, invasion, metastasis, and chemoresistance. Among all three isoforms, AKT2 is most related to cancer cell invasion, metastasis, and survival. Amplification and overexpression of AKT2 have been shown in many cancers. Accumulating evidence shows the potential role of different miRNA involvements in cancer progression by activating or suppressing AKT2 expression. The AKT2/NAB1/SPK1 pathway is a novel regulating factor of macrophage migration and cardiac remodeling after myocardial infarction. The novel mechanism of the AKT2-PKM2-STAT3/NF-kappaB axis in the

regulation of ovarian cancer progression, that both AKT2 and PKM2 may be potential targets for the treatment of ovarian cancer. AKT1 and AKT2, the AKT isoforms that are highly expressed in skeletal muscle, have distinct and overlapping functions, with AKT2 more important for insulin-stimulated glucose metabolism. In adipocytes, AKT2 versus AKT1 has greater susceptibility for insulin-mediated redistribution from cytosolic to membrane localization, and insulin also causes subcellular redistribution of AKT Substrate of 160 kDa (AS160), an AKT2 substrate and crucial mediator of insulin-stimulated glucose transport.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481