

## DGKK Protein, Human, Recombinant (His)

### General Information

Protein Construction:	Recombinant full length human DGKK was expressed by baculovirus in Sf9 insect cells using an N-terminal His tag.
Species:	Human
Expression Host:	Baculovirus-Insect Cells
Accession:	Q5KSL6
Molecular Weight:	~145 kDa

### QC Testing

Biological Activity:	The specific activity of DGKK was determined to be 26 nmol/min/mg by ADP-Glo kinase assay.
Purity:	>70% as determined by SDS-PAGE.
Formulation:	Supplied as sterile 50 mM sodium phosphate, pH 7.0, 300 mM NaCl, 150 mM imidazole, 0.1 mM PMSF, 0.25 mM DTT, 25% glycerol.

### Preparation and Storage

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

Enzymes are highly recommended to be shipped at frozen temperature with dry ice. Shipment made at ambient temperature may seriously affect the activity of the ordered products.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481