

EphA8 Protein, Human, Recombinant (GST)

General Information

| | |
|-----------------------|---|
| Protein Construction: | Recombinant human EPHA8 (565-end) was expressed by baculovirus in Sf9 insect cells using an N-terminal GST tag. |
| Species: | Human |
| Expression Host: | Baculovirus-Insect Cells |
| Accession: | P29322 |
| Molecular Weight: | ~70 kDa |

QC Testing

| | |
|----------------------|--|
| Biological Activity: | Activity has not been tested. It is theoretically active, but we cannot guarantee it. |
| Purity: | >90% as determined by SDS-PAGE. |
| Formulation: | Supplied as sterile 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 10 mM glutathione, 0.1 mM EDTA, 0.25 mM DTT, 0.1 mM PMSF, 25% glycerol. |

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Enzymes are highly recommended to be shipped at frozen temperature with dry ice. Shipment made at ambient temperature may seriously affect the activity of the ordered products.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481