

Anti-CD14 Antibody-FITC (3C10)

Product Details

Ig Type:	mouse IgG2b
Reactivity:	Human
Conjugation:	Conjugates: FITC Excitation Laser: Blue Laser (488 nm) Ex/Em: 494/520 nm
Clone:	3C10

Applications

Application:	FCM
	FCM: 5µl/Test
Recommended	Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481