

## Anti-CD14 Antibody-APC-Cy5.5 (61D3)

### Product Details

Ig Type:	mouse IgG1, $\kappa$
Reactivity:	Human
Conjugation:	Conjugates: APC-Cy5.5 Excitation Laser: Red Laser (633 nm) Ex/Em: 651/700 nm
Clone:	61D3

### Applications

Application:	FCM
Recommended	FCM: 5 $\mu$ l/Test Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x10 <sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481