

## Anti-CD14 Antibody-Biotin (61D3)

### Product Details

Ig Type: mouse IgG1,  $\kappa$   
Reactivity: Human  
Conjugation: Biotin  
Clone: 61D3

### Applications

Application: FCM  
FCM: 0.25-4 $\mu$ g/Test  
Recommended Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

### Properties

Stability & Storage: Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.  
Shipping: Shipping with blue ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481