

Anti-CD147 Antibody-PE-Cy7 (HIM6)

Product Details

| | |
|--------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Ig Type: | mouse IgG1 |
| Reactivity: | Human,Chimpanzee,Baboon,Capuchin Monkey |
| Conjugation: | Conjugates: PE-Cy7 Excitation Laser: Blue Laser (488 nm), Green Laser (532 nm), Yellow-green Laser (561 nm) Ex/Em: 565/785 nm |
| Clone: | HIM6 |

Applications

| | |
|--------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application: | FCM |
| Recommended | FCM: 5µl/Test Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test. |

Properties

| | |
|----------------------|--------------------------------------------------------------------------------|
| Stability & Storage: | Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. |
| Shipping: | Shipping with blue ice. |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481