

Anti-CD148 Antibody-PE (MEM-CD148)

Product Details

Ig Type:	mouse IgG2b
Reactivity:	Human
Conjugation:	Conjugates: PE Excitation Laser: Blue Laser (488 nm), Green Laser (532 nm), Yellow-green Laser (561 nm) Ex/Em: 566/578 nm
Clone:	MEM-CD148

Applications

Application:	FCM
Recommended	FCM: 5µl/Test Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481