

Anti-CD90 Antibody-PE (G7)

Product Details

Ig Type:	rat IgG2c, κ
Reactivity:	Mouse
Conjugation:	Conjugates: PE Excitation Laser: Blue Laser (488 nm), Green Laser (532 nm), Yellow-green Laser (561 nm) Ex/Em: 565/578 nm
Clone:	G7

Applications

Application:	FCM FCM: 5 μ l/Test
Recommended	Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481