

## Anti-CD90.2 Antibody-AF700 (30-H12)

### Product Details

Ig Type:	rat IgG2b
Reactivity:	Mouse
Conjugation:	Conjugates: AF700 Excitation Laser: Red Laser (633 nm) Ex/Em: 635/718 nm
Clone:	30-H12

### Applications

Application:	FCM
Recommended	FCM: 5µl/Test Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, $1 \times 10^6$ cells can be tested per experiment. Refer to the detailed scenario as you test.

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481