

Anti-CD94 Antibody-FITC (18d3)

Product Details

Ig Type: rat IgG2a
Reactivity: Mouse
Conjugation: Conjugates: FITC
Excitation Laser: Blue Laser (488 nm)
Ex/Em: 494/520 nm
Clone: 18d3

Applications

Application: FCM
FCM: 5µl/Test
Recommended Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Properties

Stability & Storage: Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping: Shipping with blue ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481