

Compound N026-0109

Chemical Properties


CAS No. :

Formula: C32H48O5

Molecular Weight: 512.731

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N026-0109 is a useful organic compound for research related to life sciences and the catalog number is T131308.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9503 mL	9.7517 mL	19.5034 mL
5 mM	0.3901 mL	1.9503 mL	3.9007 mL
10 mM	0.195 mL	0.9752 mL	1.9503 mL
50 mM	0.039 mL	0.195 mL	0.3901 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481