

## MBP Tag Immunomagnetic Beads

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

TargetMol MBP Tag Immunomagnetic Beads can specifically bind to proteins with a MBP tag and are suitable for immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. This product is applicable to antigen samples derived from cell lysates, cell secretion supernatants, serum, ascites, and other sources.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481