

## Isopentadecanoyl-CoA

## Chemical Properties

CAS No. :	170382-00-0
Formula:	C <sub>36</sub> H <sub>64</sub> N <sub>7</sub> O <sub>17</sub> P <sub>3</sub> S
Molecular Weight:	991.92
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Isopentadecanoyl-CoA (Isopentadecanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0081 mL	5.0407 mL	10.0815 mL
5 mM	0.2016 mL	1.0081 mL	2.0163 mL
10 mM	0.1008 mL	0.5041 mL	1.0081 mL
50 mM	0.0202 mL	0.1008 mL	0.2016 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481