

MG1 peptide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	MG1 peptide is a targeting peptide for M1 microglial cells. It can be conjugated with nanoparticles to traverse the blood-brain barrier, accurately identify M1 microglial cells, and enable controlled drug release and specific accumulation through an intelligent molecular switch. MG1 peptide is applicable in drug delivery research for ischemic stroke.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481