

P21 saRNA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	P21 saRNA is a small activating RNA (saRNA) that targets the p21 gene. It induces the expression of the p21 gene and specifically targets the p21 promoter located at the gene's transcription start site, position -322.
Targets(IC50)	CDK

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481