

LAGIPRA peptide TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

LAGIPRA peptide TFA is a long-acting GIP1R agonist. It improves insulin sensitivity by enhancing glucose processing and reducing circulating branched-chain amino acids (BCAA) and keto acids. LAGIPRA peptide TFA shows potential for research in type 2 diabetes.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481