

(Rac)-NSD2-PWWP1-IN-4

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(Rac)-NSD2-PWWP1-IN-4 (compound 7) is a potent inhibitor of NSD2-PWWP1, with IC50 and Kd values of 0.21 and 0.15 $\mu$ M, respectively. This compound is useful for tumor research.
Targets(IC50)	Histone Methyltransferase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481