

GYGGGP(GPP)5GFOGER(GPP)5GPC TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The compound GYGGGP(GPP)5GFOGER(GPP)5GPC TFA is a type I collagen mimic that can bind to integrins and be used to modify hydrogels.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481