

Volanesorsen scramble negative control

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Volanesorsenscramble negative control is the negative control form of Volanesorsen with the sequence: CAUGUTCUTCUGCATGUCAU. Volanesorsen (ISIS 304801) is an antisense oligonucleotide inhibitor targeting apolipoprotein CIII (apo-CIII) mRNA, effectively reducing triglyceride levels and improving insulin resistance.
Targets(IC50)	Apolipoprotein

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481