

Cy5-TAT acetate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5-TAT acetate is a CY5-labeled TAT peptide with excitation and emission wavelengths of 600-620 nm and 670 nm, respectively. It is suitable for use in fluorescence microscopy, particularly for single-molecule fluorescence imaging.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481