

## Human NLRP3 mRNA

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	HumanNLRP3mRNA encodes the NLRP3 protein, which functions as an inflammasome. NLRP3 is part of the NLR protein family and plays a critical role in immune responses and human diseases.
Targets(IC50)	NOD-like Receptor (NLR)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481