

DV1 TFA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DV1 TFA is a selective antagonist of the CXC chemokine receptor 4 (CXCR4). It exhibits antiviral activity by blocking the entry of HIV-1 via the CXCR4 coreceptor.
Targets(IC50)	CXCR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481