

## Sinapoyl-CoA

## Chemical Properties

CAS No. :	54429-80-0
Formula:	C32H46N7O20P3S
Molecular Weight:	973.73
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Sinapoyl-CoA (Sinapoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.027 mL	5.1349 mL	10.2698 mL
5 mM	0.2054 mL	1.027 mL	2.054 mL
10 mM	0.1027 mL	0.5135 mL	1.027 mL
50 mM	0.0205 mL	0.1027 mL	0.2054 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481