

YKYY TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	YKYY TFA is an antihypertensive peptide derived from the digest of Undaria pinnatifida (wakame), functioning as an angiotensin I-converting enzyme (ACE) inhibitor, with an IC50 of 64.2 μ M. YKYY FTA is applicable in hypertension research.
Targets(IC50)	Angiotensin-converting Enzyme (ACE)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481