

Cys-TAT-HA2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Cys-TAT-HA2 is a conjugate of Cys and TAT-HA2. It acts as a transduction compound capable of delivering proteins, such as cardiac troponin C (Cardiac troponin C), into the cytoplasm of living cells. Cys-TAT-HA2 is applicable for real-time tracking studies of exogenous proteins.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481