

Auroglaucin

Chemical Properties

CAS No. :	41451-81-4
Formula:	C19H22O3
Molecular Weight:	298.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Auroglaucin is an effective antioxidant. It exhibits activity in scavenging 1,1-diphenyl-2-picrylhydrazyl (DPPH) free radicals and superoxide radicals, with EC50 values of 74.6 and 12.3 μ M, respectively. Additionally, Auroglaucin has lower autoxidation inhibitory activity on docosahexaenoic acid (DHA).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3514 mL	16.7572 mL	33.5143 mL
5 mM	0.6703 mL	3.3514 mL	6.7029 mL
10 mM	0.3351 mL	1.6757 mL	3.3514 mL
50 mM	0.067 mL	0.3351 mL	0.6703 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481