

Vanilloyl-CoA

Chemical Properties

CAS No. :	181360-07-6
Formula:	C ₂₉ H ₄₂ N ₇ O ₁₉ P ₃ S
Molecular Weight:	917.67
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Vanilloyl-CoA (Vanilloyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0897 mL	5.4486 mL	10.8972 mL
5 mM	0.2179 mL	1.0897 mL	2.1794 mL
10 mM	0.109 mL	0.5449 mL	1.0897 mL
50 mM	0.0218 mL	0.109 mL	0.2179 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481