

STAT6 degrader-3

Chemical Properties

CAS No. :	3077464-19-5
Formula:	C47H53N11O5
Molecular Weight:	852.01
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

Biological Description

Description	STAT6degrader-3 (I-1) is a STAT6 degrader with a degradation concentration (DC50) of less than 1 nM.
Targets(IC50)	STAT

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1737 mL	5.8685 mL	11.737 mL
5 mM	0.2347 mL	1.1737 mL	2.3474 mL
10 mM	0.1174 mL	0.5868 mL	1.1737 mL
50 mM	0.0235 mL	0.1174 mL	0.2347 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481